

Product Datasheet

Human Immunoglobulin A (IgA) Antibody JP121

Product Name	Human Immunoglobulin A (IgA) Antibody
Catalog Number	JP121
Species/Host	Goat
Tested Applications	ELISA, TIA
Immunogen	Human IgA isolated from plasma
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	2-8°C
Isotype	IgG
Clonality	Polyclonal
Buffer	20 mM Sodium Phosphate, pH 7.2 with 150 mM Sodium Chloride
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.05% Sodium Azide