

## Product Datasheet

### Human Immunoglobulin A (IgA) Antibody JP121

<b>Product Name</b>	Human Immunoglobulin A (IgA) Antibody
<b>Catalog Number</b>	JP121
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, TIA
<b>Immunogen</b>	Human IgA isolated from plasma
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal, Delipidated and 0.2 µm filtered
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	20 mM Sodium Phosphate, pH 7.2 with 150 mM Sodium Chloride
<b>Specificity</b>	Single precipitin arc visible with immunoelectrophoresis
<b>Preservatives</b>	0.05% Sodium Azide