

Product Datasheet

Human Immunoglobulin M (IgM) Antibody JP126

Product Name	Human Immunoglobulin M (IgM) Antibody
Catalog Number	JP126
Species/Host	Goat
Tested Applications	TIA
Immunogen	Purified Human IgM
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	Tris-Saline, pH 7.5
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.1% Sodium Azide