

Product Datasheet

Histone H3 Antibody JP132

Product Name	Histone H3 Antibody
Catalog Number	JP132
Species/Host	Sheep
Tested Applications	IP, WB, IHC
Immunogen	Bovine histone H3-RNA complex
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term (≤ 1 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline
Cross Reactivity	Reactive with many species