

Product Datasheet

Elastase Antibody JP137

Product Name Elastase Antibody

Catalog Number JP137

Species/Host Rabbit

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Affinity Purified Polyclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Polyclonal

Purity ≥95%

Buffer Phosphate Buffered Saline, pH 7.4

Matched Pairs JP138 (capture or detection)

Preservatives 0.09% Sodium Azide