

Product Datasheet

Elastase Antibody JP137

Product Name	Elastase Antibody
Catalog Number	JP137
Species/Host	Rabbit
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	2-8°C
Clonality	Polyclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	JP138 (capture or detection)
Preservatives	0.09% Sodium Azide