

Product Datasheet

Human Lambda Light Chain Antibody JP141

Product Name	Human Lambda Light Chain Antibody
Catalog Number	JP141
Species/Host	Goat
Tested Applications	ELISA
Immunogen	Purified Human Lambda Light chain
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	TBS, pH 7.4
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.1% Sodium Azide