

## Product Datasheet

### Prealbumin Antibody JP144

<b>Product Name</b>	Prealbumin Antibody
<b>Catalog Number</b>	JP144
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Purified Human Prealbumin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal Antibody, defibrinated, delipidated, absorbed by solid phase chromatography as required and 0.2 µm filtered
<b>Storage</b>	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05M Tris, pH 7.5 with 0.5M NaCl
<b>Specificity</b>	Human
<b>Preservatives</b>	0.09% Sodium Azide