

Product Datasheet

Estradiol (E2) Antibody JP145

Product Name	Estradiol (E2) Antibody
Catalog Number	JP145
Species/Host	Rabbit
Tested Applications	ELISA, WB, IHC
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	2-8°C
Clonality	Polyclonal
Buffer	50 mM Sodium Phosphate Buffer with 0.9% NaCl, pH 7.3
Cross Reactivity	17 β Estradiol 100%; Estrone <1%; Testosterone, Corticosterone, Pregnanelone, Androstenedione, Progesterone, Estriol < 0.1%
Specificity	Estradiol
Preservatives	0.025% Proclin 300