
Product Datasheet

Transferrin Antibody JP148

Product Name	Transferrin Antibody
Catalog Number	JP148
Species/Host	Goat
Tested Applications	TIA
Immunogen	Purified Human Transferrin
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	2-8°C
Clonality	Polyclonal
Buffer	20 mM Sodium Phosphate, pH 7.2 with 150 mM Sodium Chloride
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.05% Sodium Azide