

## Product Datasheet

### Apolipoprotein CII (Apo CII) Antibody JP151

<b>Product Name</b>	Apolipoprotein CII (Apo CII) Antibody
<b>Catalog Number</b>	JP151
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Human Apo CII
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Polyclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	75 mM Sodium Phosphate, 75 mM NaCl, 0.5 mM EDTA, pH 7.2
<b>Specificity</b>	Human Apo CII
<b>Preservatives</b>	0.02% Sodium Azide