

Product Datasheet

S100B Antibody JP152

Product Name	S100B Antibody
Catalog Number	JP152
Species/Host	Rabbit
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term: 2-8°C, Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	S100b
Matched Pairs	JP152 (self-pairs)
Preservatives	0.09% Sodium Azide