

Product Datasheet

S100B Antibody JP153

Product Name S100B Antibody

Catalog Number JP153

Species/Host Rabbit

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Purified Polyclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Polyclonal

Purity >95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Specificity S100b

Matched Pairs JP153 (self-pairs)

Preservatives 0.09% Sodium Azide