

Product Datasheet

Goat Immunoglobulin G (IgG) H + L Antibody JP183

Product Name	Goat Immunoglobulin G (IgG) H + L Antibody
Catalog Number	JP183
Species/Host	Rabbit
Tested Applications	ELISA, WB, IHC
Immunogen	Purified Goat IgG
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	2-8°C
Isotype	IgG
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with non-Ig Goat serum proteins; Cross reactivity with other species IgG
Preservatives	0.05% Sodium Azide