

Product Datasheet

NT-proBNP Antibody JP190

Product Name	NT-proBNP Antibody
Catalog Number	JP190
Species/Host	Goat
Tested Applications	ELISA
Immunogen	Recombinant human NT-proBNP (full length)
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservatives	0.09% Sodium Azide