

## Product Datasheet

### Human Immunoglobulin G (IgG) Fc Fragment Antibody JP193

<b>Product Name</b>	Human Immunoglobulin G (IgG) Fc Fragment Antibody
<b>Catalog Number</b>	JP193
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	TIA
<b>Immunogen</b>	Human IgG Fc
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal, Delipidated and 0.2 µm filtered
<b>Storage</b>	2-8°C
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	200 mM Tris, pH 7.4 with 500 mM NaCl
<b>Specificity</b>	Single precipitin arc visible with immunoelectrophoresis
<b>Preservatives</b>	0.09% Sodium Azide