

## Product Datasheet

### Donkey anti-Chicken IgY (IgG) (H&L) Antibody JP194

<b>Product Name</b>	Donkey anti-Chicken IgY (IgG) (H&L) Antibody
<b>Catalog Number</b>	JP194
<b>Species/Host</b>	Donkey
<b>Immunogen</b>	Chicken IgY whole molecule
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Polyclonal Antibody
<b>Storage</b>	Short Term (6-8 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunoaffinity chromatography with agarose beads
<b>Buffer</b>	0.02M Sodium Phosphate, 0.15M Sodium Chloride, pH 7.2
<b>Preservatives</b>	0.09% Sodium Azide