

Product Datasheet

Benzoylecgonine (BE)-BSA Antigen LA004

Product Name Benzoylecgonine (BE)-BSA Antigen

Catalog Number LA004

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Liquid; Antigen-BSA conjugate

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Matched Pairs HM1207 (detection)

Preservatives 0.01% Sodium Azide