

Product Datasheet

Chloramphenicol-BSA Antigen LA005

Product Name Chloramphenicol-BSA Antigen

Catalog Number LA005

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Liquid; Antigen-BSA conjugate

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Buffer 0.1M Phosphate Buffer, pH 7.4

Preservatives 0.05% Sodium Azide