

---

## Product Datasheet

### Chloramphenicol-BSA Antigen LA005

**Product Name** Chloramphenicol-BSA Antigen

**Catalog Number** LA005

**Tested Applications** LF

**Brand** EastCoast Bio

**Form/Appearance** Liquid; Antigen-BSA conjugate

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Buffer** 0.1M Phosphate Buffer, pH 7.4

**Preservatives** 0.05% Sodium Azide