

Product Datasheet

Ethyl glucuronide (EtG)-BSA Antigen LA007

Product Name	Ethyl glucuronide (EtG)-BSA Antigen
Catalog Number	LA007
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	Distilled Water
Matched Pairs	HM128 (detection)
Preservatives	0.01% Sodium Azide