

---

## Product Datasheet

### Mephedrone-BSA Antigen LA012

**Product Name** Mephedrone-BSA Antigen

**Catalog Number** LA012

**Tested Applications** LF

**Brand** EastCoast Bio

**Form/Appearance** Liquid; Antigen-BSA conjugate

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Buffer** 10 mM Phosphate Buffered Saline

**Preservatives** None