

---

## Product Datasheet

### Morphine-BSA (higher haptens) Antigen LA016

<b>Product Name</b>	Morphine-BSA (higher haptens) Antigen
<b>Catalog Number</b>	LA016
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Buffer</b>	Distilled Water
<b>Preservatives</b>	0.05% Sodium Azide