

Product Datasheet

Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant ? 90% LA026

| | |
|----------------------------|---|
| Product Name | Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant \geq 90% |
| Catalog Number | LA026 |
| Tested Applications | ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid |
| Storage | Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Purity | \geq 90% |
| Source | |
| Buffer | 20 mM phosphate buffer, 5 mM EDTA, 0.05% SDS, pH 7.4 |
| Specificity | NS3, Core, NS4 |
| Expression System | E. coli |
| Preservatives | None |