

## Product Datasheet

### Treponema pallidum (TpN47) Antigen, Recombinant ≥85% LA029

**Product Name** Treponema pallidum (TpN47) Antigen, Recombinant ≥85%

**Catalog Number** LA029

**Tested Applications** ELISA, WB

**Alternative Names** Syphilis

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen

**Storage** -20°C. Avoid repeated freezing and thawing.

**Purity** ≥85%

**Source**

**Buffer** 25 mM Tris-HCl, pH 8.0 with 0.1% SDS, 5 mM EDTA

**Specificity** TpN47

**Expression System** E. coli