

Product Datasheet

Hepatitis C Virus (HCV-NS3) Antigen, Recombinant ≥95% LA033

Product Name Hepatitis C Virus (HCV-NS3) Antigen, Recombinant ≥95%

Catalog Number LA033

Species/Host E. coli

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage -20°C. Avoid repeated freezing and thawing.

Purity ≥90%

Source

Buffer 25 mM Tris-HCl, 0.1% SDS, 5 mM EDTA, pH 8.0

Specificity NS3

Expression System E. coli

Preservatives None