

## Product Datasheet

### Hepatitis C Virus (HCV-NS3) Antigen, Recombinant ≥95% LA033

**Product Name** Hepatitis C Virus (HCV-NS3) Antigen, Recombinant ≥95%

**Catalog Number** LA033

**Species/Host** E. coli

**Tested Applications** ELISA

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (His-tag)

**Storage** -20°C. Avoid repeated freezing and thawing.

**Purity** ≥90%

**Source**

**Buffer** 25 mM Tris-HCl, 0.1% SDS, 5 mM EDTA, pH 8.0

**Specificity** NS3

**Expression System** E. coli

**Preservatives** None