

Product Datasheet

Hepatitis C Virus (HCV-core) Antigen, Recombinant >98% LA034

Product Name	Hepatitis C Virus (HCV-core) Antigen, Recombinant >98%
Catalog Number	LA034
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	$\geq 95\%$
Source	
Buffer	25 mM Tris-HCl, 0.1% SDS, 5 mM EDTA, pH 8.0
Specificity	Core
Expression System	E. coli
Preservatives	None