

Product Datasheet

Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant ? 95% LA035

Product Name	Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant \geq 95%
Catalog Number	LA035
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	\geq 95 %
Source	
Buffer	25 mM Tris-HCl, 500 mM NaCl, 3 mM EDTA, pH 7.5
Specificity	NS3, Core, NS4, NS5
Matched Pairs	LA036 (detection)
Expression System	E. coli