

## Product Datasheet

### Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant ? 95% LA035

<b>Product Name</b>	Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant $\geq$ 95%
<b>Catalog Number</b>	LA035
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	$\geq$ 95 %
<b>Source</b>	
<b>Buffer</b>	25 mM Tris-HCl, 500 mM NaCl, 3 mM EDTA, pH 7.5
<b>Specificity</b>	NS3, Core, NS4, NS5
<b>Matched Pairs</b>	LA036 (detection)
<b>Expression System</b>	E. coli