

## Product Datasheet

### Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant >90% LA036

**Product Name** Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant >90%

**Catalog Number** LA036

**Tested Applications** LF

**Brand** EastCoast Bio

**Form/Appearance** Liquid

**Storage** -20°C

**Purity** >90%

**Source**

**Buffer** 20 mM Tris-HCl, 160 mM NaCl, 3 mM EDTA, pH 8.0 with 1 mM DTT, and 20% glycerol

**Specificity** NS3, Core, NS4

**Matched Pairs** LA035 (capture)

**Expression System** E. coli