

Product Datasheet

Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant >90% LA036

Product Name Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant >90%

Catalog Number LA036

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Liquid

Storage -20°C

Purity >90%

Source

Buffer 20 mM Tris-HCl, 160 mM NaCl, 3 mM EDTA, pH 8.0 with 1 mM DTT, and 20% glycerol

Specificity NS3, Core, NS4

Matched Pairs LA035 (capture)

Expression System E. coli