

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1) Antigen, Recombinant LA038

Product Name	Human Immunodeficiency Virus 1 (HIV 1) Antigen, Recombinant
Catalog Number	LA038
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	
Buffer	20 mM Phosphate Buffer, 0.1% SDS, pH 7.4
Specificity	gp120 and gp41
Molecular Weight	43 kDa
Matched Pairs	LA249, LA250 or LA779 (capture)
Expression System	E. coli