
Product Datasheet

Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >90% LA039

Product Name	Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >90%
Catalog Number	LA039
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	2-8°C
Source	
Buffer	20 mM Phosphate Buffer, 0.1% SDS, pH 7.4
Specificity	gp36
Matched Pairs	LA040
Expression System	E. coli