

Product Datasheet

Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >98% LA040

Product Name	Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >98%
Catalog Number	LA040
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	2-8°C
Purity	>98%
Source	
Buffer	20 mM Phosphate Buffer, pH 7.4
Specificity	gp36
Matched Pairs	LA039
Expression System	E. coli