

Product Datasheet

Treponema pallidum (T. pallidum) Chimeric Antigen, Recombinant >90% LA041

Product Name	Treponema pallidum (T. pallidum) Chimeric Antigen, Recombinant >90%
Catalog Number	LA041
Tested Applications	LF
Alternative Names	Syphilis
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	≥90%
Source	
Buffer	25 mM Tris-HCl, pH 8.3 with 0.1% SDS
Specificity	TpN15, TpN17, TpN47
Matched Pairs	LA042 (capture)
Expression System	E. coli