

Product Datasheet

Dengue Virus (DENV-Env Serotype 3) Antigen, Recombinant ≥90% LA046

Product Name	Dengue Virus (DENV-Env Serotype 3) Antigen, Recombinant ≥90%
Catalog Number	LA046
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen (His-tag)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	≥90%
Buffer	20 mM Phosphate Buffer, 50 mM NaCl, 0.1% ProClin 300, pH 7.4
Specificity	Dengue III, Envelope
Expression System	E. coli