

Product Datasheet

Dengue Virus (DENV-Env Serotype 1) Chimeric Antigen, Recombinant ≥90% LA047

Product Name	Dengue Virus (DENV-Env Serotype 1) Chimeric Antigen, Recombinant ≥90%
Catalog Number	LA047
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	2-8°C
Purity	≥90%
Buffer	20 mM phosphate buffer, pH 7.4 with 0.15M NaCl
Specificity	Dengue I, Envelope (epitope 250-376 a.a.)
Expression System	E. coli
Preservatives	0.1% Proclin 300