

Product Datasheet

Troponin I , Cardiac Antigen, Recombinant ≥98% LA146

Product Name	Troponin I , Cardiac Antigen, Recombinant ≥98%
Catalog Number	LA146
Brand	EastCoast Bio
Form/Appearance	Liquid, Recombinant Antigen
Storage	Short Term (< 1 week): -20°C. Long Term: ≤-70°C. Avoid repeated freezing and thawing.
Purity	≥98%
Buffer	6M Urea, 50mM Tris-HCl, pH 8.0
Expression System	E. coli