

## Product Datasheet

### Human Immunoglobulin E (IgE) Antigen, Native ≥95% LA170

<b>Product Name</b>	Human Immunoglobulin E (IgE) Antigen, Native ≥95%
<b>Catalog Number</b>	LA170
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	≤-20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Testing: Shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.
<b>Purity</b>	≥95% by SDS-PAGE
<b>Source</b>	Human Myeloma Plasma
<b>Buffer</b>	15 mM potassium phosphate, pH 7.4, 150 mM NaCl and 0.05% NaN <sub>3</sub>