

Product Datasheet

Human Immunoglobulin G (IgG), Fc Fragment, Native ≥95% LA172

Product Name	Human Immunoglobulin G (IgG), Fc Fragment, Native ≥95%
Catalog Number	LA172
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-20°C. Avoid repeated freezing and thawing.
Note	Testing: Shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.
Purity	≥95% by SDS-PAGE
Source	Human Plasma
Buffer	50 mM Tris-HCl, pH 8.0 with 200 mM NaCl
Preservatives	0.05% Sodium Azide