

## Product Datasheet

### Ethyl glucuronide (EtG)-BSA Antigen LA189

<b>Product Name</b>	Ethyl glucuronide (EtG)-BSA Antigen
<b>Catalog Number</b>	LA189
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid Antigen Conjugate
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Buffer</b>	10 mM Sodium Phosphate, pH 7.4 with 0.15M NaCl
<b>Matched Pairs</b>	Anti-EtG R12403 SPTNZ-1 (detection)
<b>Preservatives</b>	None