

Product Datasheet

Ethyl glucuronide (EtG)-BSA Antigen LA189

| | |
|----------------------------|---|
| Product Name | Ethyl glucuronide (EtG)-BSA Antigen |
| Catalog Number | LA189 |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid Antigen Conjugate |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Buffer | 10 mM Sodium Phosphate, pH 7.4 with 0.15M NaCl |
| Matched Pairs | Anti-EtG R12403 SPTNZ-1 (detection) |
| Preservatives | None |