

Product Datasheet

Legionella pneumophila Antigen, Native LA205

| | |
|----------------------------|--|
| Product Name | Legionella pneumophila Antigen, Native |
| Catalog Number | LA205 |
| Tested Applications | ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid |
| Storage | ≤-70°C |
| Note | Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious. |
| Source | Strain ATCC 33153 |
| Buffer | 20 mM Phosphate Buffer, pH 7.4 |
| Preservatives | None |