

## Product Datasheet

### Cytomegalovirus (CMV) Antigen, Native LA209

<b>Product Name</b>	Cytomegalovirus (CMV) Antigen, Native
<b>Catalog Number</b>	LA209
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	≤-70°C
<b>Note</b>	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
<b>Source</b>	MRC-5 Cells/Strain AD169
<b>Buffer</b>	10 mM Phosphate, pH 7.4 with 140 mM NaCl and <0.1% CHAPS