

Product Datasheet

C-Reactive Protein, Canine (cCRP) Antigen, Native ≥98% LA216

Product Name	C-Reactive Protein, Canine (cCRP) Antigen, Native ≥98%
Catalog Number	LA216
Tested Applications	LF, ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Liquid, affinity purified
Storage	≤-20°C. Avoid repeated freezing and thawing.
Purity	≥98% (20 kDa bands with SDS-PAGE)
Source	Dog serum
Buffer	100 mM Tris HCl, pH 7.5 with 200 mM NaCl, and 2 mM CaCl ₂
Specificity	Strongly reactive with tested Canine CRP antibodies
Preservatives	0.05% Sodium Azide