

Product Datasheet

Calicivirus, Feline Antigen, Native LA223

| | |
|----------------------------|--|
| Product Name | Calicivirus, Feline Antigen, Native |
| Catalog Number | LA223 |
| Tested Applications | ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid; Native Antigen |
| Storage | ≤ 70°C |
| Note | Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious. |
| Source | FC 3TG Cells |
| Buffer | 50 mM Glycine, pH 9 with 140 mM NaCl and <0.1% CHAPS |
| Preservatives | None |