

Product Datasheet

Adenovirus Type 2, Canine Native Antigen LA224

Product Name	Adenovirus Type 2, Canine Native Antigen
Catalog Number	LA224
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Native Antigen
Storage	≤ -70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	A-72 Cells/Virus Serotype 2
Buffer	50 mM Glycine, pH 9 with 140 mM NaCl and <0.1% CHAPS
Preservatives	None