

## Product Datasheet

### Adenovirus, Mouse Native Antigen LA230

<b>Product Name</b>	Adenovirus, Mouse Native Antigen
<b>Catalog Number</b>	LA230
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Native Antigen
<b>Storage</b>	≤ -70°C
<b>Note</b>	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
<b>Source</b>	NIH 3T3 Cells/Virus Strain FL
<b>Buffer</b>	50 mM Glycine, pH 9 with 140 mM NaCl, <0.3% CHAPS and <0.01% SDS
<b>Preservatives</b>	None