

## Product Datasheet

### Pneumonia Virus of Mice (PVM) Antigen, Native LA231

<b>Product Name</b>	Pneumonia Virus of Mice (PVM) Antigen, Native
<b>Catalog Number</b>	LA231
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	≤-70°C
<b>Note</b>	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
<b>Source</b>	BHK-21 cells/Strain No. 15
<b>Buffer</b>	50 mM Glycine, pH 9 with <0.2% CHAPS and 0.01% SDS
<b>Preservatives</b>	None